# CHEMISTRY 235 EXPERIMENT 5 QUALITATIVE ANALYSIS

## METHOD AND RESULTS

The compounds A, B, and C contain the ions  $Co^{2+}$ ,  $Ba^{2+}$  and  $Ca^{2+}$ . Identify the ion present in each compound. Write the equation, where relevant, for each reaction in the deduction column.

#### **COMPOUND A**

TEST	OBSERVATIONS	DEDUCTIONS
Make an aqueous soln. of A for the following tests. Use a fresh portion for each test unless otherwise instructed.		
a) Add NaOH(aq) until in excess.		
b) Add NH <sub>3</sub> (aq) until in excess.		
c) Add Na <sub>2</sub> CO <sub>3</sub> (aq)		
d) Add dil. HCl then $H_2S(aq)$ .		
e) Add ethanol followed by NH <sub>4</sub> SCN(aq).		

#### **COMPOUND B**

TEST	OBSERVATIONS	DEDUCTIONS
a) Perform a flame test on B.		
Prepare a solution of B and use it for the following tests. Use a fresh portion for each test unless otherwise instructed.		

TEST	OBSERVATIONS	DEDUCTIONS
b) Add NaOH(aq) until in excess.		
c) Add NH <sub>3</sub> (aq) until in excess.		
d) Add Na <sub>2</sub> CO <sub>3</sub> (aq) and then dil. HCl or dil. HNO <sub>3</sub> .		
e) Add K <sub>2</sub> SO <sub>4</sub> (aq) then dil. HCl or dil. HNO <sub>3</sub> .		
f) Add Na <sub>3</sub> PO <sub>4</sub> (aq) then dil. HNO <sub>3</sub> .		
g) Add Na <sub>2</sub> (COO) <sub>2</sub> (aq) <sup>1</sup> then dil. HCl or dil. HNO <sub>3</sub> .		
h) Add H <sub>2</sub> S soln.		
i) Add K <sub>2</sub> CrO <sub>4</sub> (aq) and divide result into 2 parts.		
<ul><li>α) To one part</li><li>add dil HCl or dil.</li><li>HNO<sub>3</sub>.</li></ul>		
$\beta$ ) To 2nd part add ethanoic acid.		

<sup>&</sup>lt;sup>1</sup>Sodium oxalate (sodium ethanedioate) solution

### COMPOUND C

TESTS	OBSERVATIONS	DEDUCTIONS
a) Perform a flame test on C.		
Prepare an aqueous soln. of C and use it for the following tests.		
b) Add NaOH(aq) until in excess.		
c) Add NH <sub>3</sub> (aq) until in excess.		
d) Add Na <sub>2</sub> CO <sub>3</sub> (aq) then dil. HCl or dil. HNO <sub>3</sub> .		
e) Add Na <sub>3</sub> PO <sub>4</sub> (aq) then dil. HCl or dil. HNO <sub>3</sub> .		
f) Add Na <sub>2</sub> (COO) <sub>2</sub> (aq) then dil. HCl or dil. HNO <sub>3</sub> .		
g) Add $K_2CrO_4(aq).$		